

Tests available

Pathogen	Sample collection
<ul style="list-style-type: none">• <i>Borrelia</i> spp.	Blood
<ul style="list-style-type: none">• <i>Theileria equi</i>	
<ul style="list-style-type: none">• <i>Babesia caballi</i>	
<ul style="list-style-type: none">• <i>Anaplasma phagocytophylum</i>	
<ul style="list-style-type: none">• <i>Leptospira</i> pathogenic species	
<ul style="list-style-type: none">• EHV-1	Nasopharyngeal swab
<ul style="list-style-type: none">• EHV-4	
<ul style="list-style-type: none">• EHV-1 and EHV-4	
<ul style="list-style-type: none">• Equine influenza	
<ul style="list-style-type: none">• Respiratory panel: EHV-1, EHV-4, equine influenza, <i>Str. equi equi</i>	
<ul style="list-style-type: none">• Respiratory-virus panel: EHV-1 and EHV-4, equine influenza	
<ul style="list-style-type: none">• <i>Str. equi equi</i>. Has greater sensitivity for this organism than the respiratory-panel test	